



FIND THEIR PERSONALIZED PATH WITH **LIQUID BIOPSY**

See a way forward for all patients with solid tumours – using FoundationOne®Liquid CDx, a **liquid biopsy** comprehensive genomic profiling service, as a **minimally-invasive** option, **alternative or complementary** to FoundationOne®CDx, at **optimal times beneficial to their treatment journey**¹⁻³

© 2020 Foundation Medicine, Inc., Foundation Medicine® and FoundationOne® are trade-marks of Foundation Medicine, Inc. Content modified by Hoffmann-La Roche Limited for Canadian use.


See more, do more



FOUNDATION
MEDICINE®



FoundationOne Liquid CDx provides insights that inform research or treatment decisions for individual patients across all solid tumours^{2,3}



*309 genes with complete exon coverage, 15 genes with select intronic or non-coding regions only.


**FoundationOne Liquid CDx reports MSI-H status.

FoundationOne Liquid CDx has been validated with >7,500 samples covering >30,000 unique variants across 324 genes and 37 cancer indications^{†‡}


- 150 X more samples than the typical number of patient samples for laboratory test validation^{†6} and similar to the number of samples used for validation of FoundationOne CDx⁷
- 12 X more variants than the number of variants used for validation of our previous liquid-based test FoundationOne Liquid⁸

[†]Validation against an externally validated cfDNA assay.

[‡]According to New York State guidance.



FoundationOne Liquid CDx reports Tumour Fraction, the percentage of tumour-derived DNA in the circulating cell-free DNA, providing you with a measure of the likelihood of detecting genomic alterations in the liquid biopsy sample^{2,3}



- FoundationOne Liquid CDx reports Tumour Fraction estimates above ~5%^{2,3,9}
- The Tumour Fraction measurement is designed as a decision insight that can guide reflex testing decisions
- A higher Tumour Fraction indicates a higher concordance to tissue-based testing^{2,9}

bTMB, blood Tumour Mutational Burden; CDx, companion diagnostic; cfDNA, circulating cell-free DNA; MSI, Microsatellite Instability; TF, Tumour Fraction.

Quick and convenient single blood draw helps avoid invasive biopsies and enables faster treatment decisions^{2,3,10-15}




Two tubes of whole blood (8.5 mL each)




Turnaround time from receipt of sample at our laboratory to report

Typical timelines associated with liquid- versus tissue-based CGP














- Liquid biopsy genomic profiling may help reduce the distress which patients may experience from a longer waiting time for scheduling a biopsy and receiving the results of tissue-based testing^{11,12}
- FoundationOne Liquid CDx may also help reduce requirements for healthcare infrastructure and resources compared with tissue-based testing¹⁵
- Phlebotomy services for FoundationOne Liquid CDx are available through Bayshore Infusion Clinics. Please email canada.fmi@roche.com for more information.

Suitable for all solid tumours based on your patient's profile and the available clinical evidence^{2,3,13,14}



Part of our proven portfolio of high-quality, comprehensive genomic profiling services for all solid tumours, supporting current and future personalized clinical and research decisions at optimal times beneficial to your patients' treatment journeys^{1-3,16,17}

Primary diagnosis and therapy initiation	Monitoring	Progression and further lines of therapy selection based on new/evolving resistance mutations
<div><div> or </div><div>TISSUE (if available) LIQUID</div><div>Use either FoundationOne CDx or FoundationOne Liquid CDx based on medical need and available evidence; reflex to the other service if sample is inadequate or no actionable genomic alterations detected</div></div> <div><div> + </div><div>Complement FoundationOne CDx with FoundationOne Liquid CDx to capture the full disease heterogeneity</div></div> <div>Tissue-based testing is often recommended at diagnosis due to availability of tissue</div>	<div><div>  </div><div>LIQUID</div><div>FoundationOne Liquid CDx may also be used to track responses to therapy and for early diagnosis of relapse</div></div>	<div><div> or </div><div>LIQUID TISSUE (if available)</div><div>Use FoundationOne Liquid CDx and reflex to FoundationOne CDx if no actionable genomic alterations detected in the blood sample</div></div> <div><div> + </div><div>Complement FoundationOne Liquid CDx with FoundationOne CDx to capture the full disease heterogeneity</div></div> <div>Liquid-based testing is often suitable at progression due to frequent challenges with tissue biopsy and the increase in ctDNA as the disease progresses</div>

*FoundationOne Liquid CDx reports MSI-H status.

CDx, companion diagnostic; ctDNA: circulating tumour DNA.



FoundationOne Liquid CDx sequences cfDNA which can originate from primary and metastatic tumour sites, thereby capturing the full disease heterogeneity across the body as the disease evolves^{2,3,13,14}

Tumour tissue


Apoptosis
Necrosis secretion

Healthy tissue

Main source of cell-free DNA



Tumour heterogeneity within the primary tumour



Tumour heterogeneity across the body

FoundationOne Liquid CDx can help inform therapy selection at progression^{13,14,18}

- Tumour heterogeneity is associated with the development of new and evolving resistance mutations and relapse¹⁸
- FoundationOne Liquid CDx is often suited at progression owing to the increase in ctDNA levels and clonal evolution of the tumour as the disease progresses^{13,14}




Figure adapted from reference 18.

FoundationOne Liquid CDx is suited for NSCLC, breast cancer, ovarian cancer, prostate cancer, and other solid tumours^{2,3,19-39}

- FoundationOne Liquid CDx covers professional guidelines, and other clinically relevant genomic alterations and signatures* for NSCLC, breast, ovarian and prostate cancer^{2,3,19-31,33-39}



NSCLC



BREAST





PROSTATE



OVARIAN

*FoundationOne Liquid CDx reports MSI-H status.
CDx, companion diagnostic; ctDNA: circulating tumour DNA; NSCLC, non-small cell lung cancer.

In-depth report helps inform clinical decision-making⁴⁰



- Consistency in the reports across different FoundationOne services aids comparison of the results across the patient journey^{40,41}

*FoundationOne Liquid CDx reports MSI-H status

†For additional information on the NCCN categories please refer to the NCCN Compendium* (www.nccn.org).

bTMB, blood Tumour Mutational Burden; CDx, companion diagnostic; MSI, Microsatellite Instability; NCCN, National Comprehensive Cancer Network; VAF, variant allele frequency

FoundationOne Liquid CDx is able to detect novel and known variants of all four main classes of genomic alterations in 324 cancer-related genes in all solid tumours^{2,3}

Genes baited with enhanced sensitivity[†] (selected based on increased actionability with current or future targeted therapies) are highlighted^{2,3}

ABL1 [Exons 4–9]	CALR	CYP17A1	FGFR4	KDM6A	MYCL (MYCL1)	POLD1	SMAD4
ACVR1B	CARD11	DAXX	FH	KDR	MYCN	POLE	SMARCA4
AKT1 [Exon 3]	CASP8	DDR1	FLCN	KEAP1	MYD88 [Exon 4]	PPARG	SMARCB1
AKT2	CBFB	DDR2 [Exons 5,17,18]	FLT1	KEL	NBN	PPP2R1A	SMO
AKT3	CBL	DIS3	FLT3 [Exons 14,15,20]	KIT [Exons 8,9,11,12,13,17, Intron 16]	NF1	PPP2R2A	SNCAIP
ALK [Exons 20–29, Introns 18,19]	CCND1	DNMT3A	FOXL2	KLHL6	NF2	PRDM1	SOC31
ALOX12B	CCND2	DOT1L	FUBP1	KMT2A (MLL) [Introns 6,8–11, Intron 7]	NFE2L2	PRKARIA	SOX2
AMER1 (FAM123B)	CCND3	EED	GABRA6	KMT2D (MLL2)	NFKBIA	PRKCI	SOX9
APC	CCNE1	EGFR [Introns 7,15,24–27]	GATA3	KRAS	NKX2–1	PTCH1	SPEN
AR	CD22	EP300	GATA4	LTK	NOTCH1	PTEN	SPOP
ARAF [Exons 4,5,7, 11,13,15,16]	CD70	EPHA3	GATA6	LYN	NOTCH2 [Intron 26]	PTPN11	SRC
ARFRP1	CD79A	EPHB1	GNA11 [Exons 4,5]	MAF	NOTCH3	PTPRO	STAG2
ARID1A	CD79B	EPHB4	GNA13	MAP2K1 (MEK1) [Exons 2,3]	NPM1 [Exons 4–6,8,10]	QKI	STAT3
ASXL1	CD274 (PD-L1)	ERBB2	GNAQ [Exons 4,5]	MAP2K2 (MEK2) [Exons 2–4,6,7]	NRAS [Exons 2,3]	RAC1	STK11
ATM	CDC73	ERBB3 [Exons 3,6,7, 8,10,12,20, 21,23,24,25]	GNAS [Exons 1,8]	MAP2K4	NSD3 (WHSC1L1)	RAD21	SUFU
ATR	CDH1	ERBB4	GRM3	MAP3K1	NTSC2	RAD51	SYK
ATRX	CDK12	ERCC4	GSK3B	MAP3K13	NTRK1 [Exons 14,15, Introns 8–11]	RAD51B	TBX3
AURKA	CDK4	ERG	H3F3A	MAPK1	NTRK2 [Intron 12]	RAD51C	TEK
AURKB	CDK6	ERRF1	HDAC1	MCL1	NTRK3 [Exons 16,17]	RAD51D	TERC {ncRNA}
AXIN1	CDK8	ESR1 [Exons 4–8]	HGF	MDM2	P2RY8	RAD52	TERT {Promoter}
AXL	CDKN1A	EZH2 [Exons 4,16,17,18]	HNFA	MDM4	PALB2	RAD54L	TET2
BAP1	CDKN1B	FAM46C	HRAS [Exons 2,3]	MED12	PARK2	RAF1 [Exons 3,4,6,7,10,14,15,17, Introns 4–8]	TGFBR2
BARD1	CDKN2A	FANCA	HSD3B1	MEF2B	PARP1	RARA [Intron 2]	TIPARP
BCL2	CDKN2B	FANCC	ID3	MEN1	PARP2	RB1	TNFAIP3
BCL2L1	CDKN2C	FANCG	IDH1 [Exon 4]	MERTK	PARP3	RBM10	TNFRSF14
BCL2L2	CEBPA	FANCL	IDH2 [Exon 4]	MET	PAX5	REL	TP53
BCL6	CHEK1	FAS	IGF1R	MITF	PBRM1	RET [Introns 7,8, Exons 11,13–16, Introns 9–11]	TSC1
BCOR	CHEK2	FBXW7	IKBKE	MKNK1	PDCD1 (PD-1)	RICTOR	TSC2
BCORL1	CIC	FGF10	IKZF1	MLH1	PDCD1LG2 (PD-L2)	RNF43	TYRO3
BRAF [Exons 11–18, Introns 7–10]	CREBBP	FGF12	INPP4B	MPL [Exon 10]	PDGFRA [Exons 12,18, Introns 7,9,11]	ROS1 [Exons 31,36–38,40, Introns 31–35]	U2AF1
BRCA1 [Introns 2,7, 8,12,16,19,20]	CRKL	FGF14	IRF2	MRE11A	PDGFRB [Exons 12–21,23]	RPTOR	VEGFA
BRCA2 [Intron 2]	CSF1R	FGF19	IRF4	MSH2 [Intron 5]	PDK1	SDHA	VHL
BRD4	CSF3R	FGF23	IRS2	MSH3	PIK3C2B	SDHB	WHSC1
BRIPI	CTCF	FGF3	JAK1	MSH6	PIK3C2G	SDHC	WT1
BTG1	CTNNA1	FGF4	JAK2 [Exons 14]	MST1R	PIK3CA [Exons 2,3,5–8,10,14,19,21] (Coding Exons 1,2, 4–7,9,13,18,20)	SDHD	XPO1
BTG2	CTNNB1 [Exon 3]	FGF6	JAK3 [Exons 5,11,12,13,15,16]	MTAP	PIK3CB	SETD2	XRCC2
BTK [Exons 2,15]	CUL3	FGFR1 [Introns 1,5, Intron 17]	JUN	MTOR [Exons 19,30,39,40,43–45,47,48,53,56]	PIK3R1	SF3B1	ZNF217
C11orf30 (EMSY)	CUL4A	FGFR2 [Intron 1, Intron 17]	KDM5A	MUTYH	PIM1	SGK1	ZNF703
C17orf39 (GID4)	CXCR4	FGFR3 [Exons 7,9 (alternative designation exon 10),14,18, Intron 17]	KDM5C	MYC [Intron 1]	PMS2	SMAD2	

*Overall 309 genes with complete coding exonic coverage, 15 genes with selected intronic or non-coding regions only.

[†]**75 genes are baited for enhanced sensitivity** using ultra-deep sequencing coverage; other genomic regions are baited for high sensitivity using high sequencing coverage.

[†]Genes with selected non-coding coverage only; additional genes with selected non-coding coverage only are: **BCR** (introns 8,13,14); **CD74** (introns 8–6); **ETV4** (introns 5,6); **ETV5** (introns 6,7); **ETV6** (introns 5,6); **EWSR1** (introns 7,13); **EZR** (introns 9–11); **MYB** (intron 14); **NUTM1** (intron 1); **RSPO2** (intron 1); **SDC4** (intron 2); **SLC34A2** (intron 4); **TMPPRSS2** (introns 1–3).

CDx, companion diagnostic; ncRNA, non-coding RNA.

Use FoundationOne Liquid CDx, our liquid biopsy comprehensive genomic profiling service, as a minimally invasive option, alternative or complementary to FoundationOne CDx, for all patients with solid tumours at optimal times beneficial to their treatment journey¹⁻³



Comprehensive panel

- Assesses all four main classes of genomic alterations in 324 genes plus MSI* and bTMB
- Complete exon coverage of 309 genes
- Additional coverage of selected intronic or non-coding regions
- Enhanced sensitivity for 75 genes associated with increased actionability with targeted therapies



Tumour Fraction

- Reports Tumour Fraction as a measure of the likelihood of detecting genomic alterations in the liquid biopsy sample



Extensively analytically and clinically validated comprehensive platform

- High sensitivity and specificity[†] for key genomic alterations, MSI* and bTMB across all solid tumours
- Analytical and clinical validation based on >7,500 samples covering >30,000 unique variants in 324 genes and 37 cancer indications[‡]



Faster treatment decisions

- Quick and convenient single blood draw and improved turnaround time of <2 weeks to enable faster treatment decisions
- Helps avoid invasive biopsies



Supporting personalized decision-making across the patient journey

- Part of our proven portfolio of high-quality, comprehensive genomic profiling services for all solid tumours, supporting current and future personalized clinical and research decisions at optimal times beneficial to your patients' treatment journeys
- Option of reflex testing to FoundationOne CDx if no actionable genomic alterations are detected in the liquid biopsy sample



High accuracy

- Accurately detects the main types of genomic alterations, MSI* and bTMB at very low allele frequencies, Tumour Fractions or percent unstable loci
- Low overall false positive rate of 1 in 8000

*FoundationOne Liquid CDx reports MSI-H status.

[†]75 genes are baited with enhanced sensitivity for all variant types (selected based on increased actionability with current or future targeted therapies; for more information of these 75 genes, please refer to our full gene list); other genomic regions are baited with high sensitivity.

[‡]Validation against an externally validated cfDNA assay.

bTMB: blood Tumour Mutational Burden; CDx, companion diagnostic; MSI, Microsatellite Instability; TAT, turnaround time; TF, Tumour Fraction.

Contacts

For more information, please contact Roche at canada.fmi@roche.com

References

1. FoundationOne Liquid CDx FDA Approval, 2020. Available at: <https://www.foundationmedicine.com/press-releases/445c1f9e-6cbb-488b-84ad-5f133612b721> (Accessed August 2020). 2. Data on file: FoundationOne Liquid CDx Technical Specifications, 2020. Available at: www.eifu.online/FMI/190070862 (Accessed August 2020). 3. Data on file: Clinical and analytical validation data file for FoundationOne Liquid CDx. 4. FDA approves pembrolizumab for adults and children with TMB-H solid tumors, 2020. Available at: <https://www.fda.gov/drugs/drug-approvals-and-databases/fda-approves-pembrolizumab-adults-and-children-tmb-h-solid-tumors> (Accessed August 2020). 5. Marabelle A et al. *Ann Oncol*. 2019;30 (suppl. 5):v475-v532. 6. Next Generation Sequencing (NGS) Guidelines For Somatic Genetic Variant Detection, 2018. Available at: <https://www.wadsworth.org/sites/default/files/WebDoc/3NextGenSeqONCOGuidelines%2012318.pdf> (Accessed August 2020). 7. FoundationOne®CDx Clinical Validation, 2017. Available at: <http://www.foundationmedicine.com/genomic-testing/foundation-one-cdx> (Accessed August 2020). 8. Clark TA et al. *J Mol Diagn* 2018; 20: 686-702. 9. Foundation Medicine; data on file 10. Data on file: FoundationOne Liquid CDx Specimen instructions. 11. Lebel S et al. *J Psychosom Res* 2003; 55: 437-443. 12. Hayes Balmadrid MA et al. *J Health Psychol* 2017; 22:561-571. 13. Francis G, Stein S. *Int J Mol Sci* 2015; 16: 14122-14142. 14. Siravegna G et al. *Ann Oncol* 2019; 30: 1580-1590. 15. Arnaud A. *Value Health* 2016; 19: PCN57. 16. FoundationOne CDx Technical Specifications, 2018. Available at: www.roche.foundationmedicine.com/ficdxtech (Accessed August 2020). 17. FoundationOne CDx FDA Approval, 2017. Available at: <https://www.fda.gov/medicaldevices/recently-approved-devices/foundationone-cdx-p170019> (Accessed August 2020). 18. Lim ZF, Ma PC. *J Hematol Oncol* 2019; 12: 134. 19. Planchar D et al. *Ann Oncol* 2019; 30: 863-870. 20. NCCN Clinical Practice Guidelines in Oncology. Non-Small Cell Lung Cancer. Version 6.2020, June 2020. Available at: https://www.nccn.org/professionals/physician_gls/pdf/nscl.pdf (Accessed August 2020). 21. Cardoso F et al. *Ann Oncol* 2019; 30: 1674. 22. NCCN Clinical Practice Guidelines in Oncology. Breast Cancer. Version 5.2020, July 2020. Available at: https://www.nccn.org/professionals/physician_gls/pdf/breast.pdf (Accessed August 2020). 23. NCCN Clinical Practice Guidelines in Oncology. Prostate Cancer. Version 2.2020, May 2020. Available at: https://www.nccn.org/professionals/physician_gls/pdf/prostate.pdf (Accessed August 2020). 24. Parker C et al. *Ann Oncol* 2020; 31: 1119-1134. 25. NCCN Clinical Practice Guidelines in Oncology. Ovarian Cancer. Version 1.2020, March 2020. https://www.nccn.org/professionals/physician_gls/pdf/ovarian.pdf (Accessed August 2020). 26. Colombo N et al. *Ann Oncol* 2019; 30: 672-705. 27. Mateo J et al. *Eur Urol* 2017; 71: 417-425. 28. Khagi Y et al. *Clin Cancer Res* 2017; 23: 5729-5736. 29. Socinski M. *Ann Oncol* 2019; 30(suppl. 5): v851-v934. 30. Gandara DR et al. *Nat Med* 2018; 24: 1441-1448. 31. Abida W et al. *JAMA Oncol* 2019; 5: 471-478. 32. Conway AM et al. *Br J Cancer* 2019; 120: 141-153. 33. Amatu A et al. *ESMO Open* 2016; 1: e000023. 34. Kok M et al. *ESMO Open* 2019; 4(Suppl 2): e000511. 35. Zhao P et al. *J Hematol Oncol* 2019; 12: 54. 36. Pikor LA et al. *Lung Cancer* 2013; 82:179-189. 37. Tsao AS et al. *J Thorac Oncol* 2016; 11: 613-638. 38. Sanchez-Vega F et al. *Cell* 2018; 173: 321-37.e10. 39. Kroeger PT Jr, et al. *Curr Opin Obstet Gynecol* 2017; 29: 26-34. 40. Data on file; FoundationOne Liquid CDx mock report, 2020. 41. FoundationOne CDx Sample Report. Available at: [https://www.foundationmedicine.nl/content/dam/rfm/sample-reports/ficdx/eu_version_-_ema_without_page_1/FICDX%20EU%20Sample%20Report%20\(Lung\).pdf](https://www.foundationmedicine.nl/content/dam/rfm/sample-reports/ficdx/eu_version_-_ema_without_page_1/FICDX%20EU%20Sample%20Report%20(Lung).pdf) (Accessed August 2020).

If you require this information in an accessible format, please contact Roche at 1-800-561-1759.

© 2020 Foundation Medicine, Inc., Foundation Medicine® and FoundationOne® are registered trade-marks of Foundation Medicine, Inc., used under license. Content modified by Hoffmann-La Roche Limited for Canadian use.

FoundationOne®Liquid CDx and FoundationOne®CDx is a next-generation sequencing service performed by Foundation Medicine, Inc. at laboratories located in Cambridge, Massachusetts, USA. Hoffmann-La Roche Limited markets this service in Canada under an exclusive license from Foundation Medicine, Inc. The FoundationOne®Liquid CDx reports are prepared exclusively by Foundation Medicine, Inc., without Hoffmann-La Roche Limited's involvement. The reports may include scientific information about experimental medicinal product(s) or uses that are not approved and/or available in Canada. Hoffmann-La Roche Limited does not endorse or recommend the unapproved use of any medicinal product.



**FOUNDATION
MEDICINE®**

